



MEDICA PONTEM

DIADEZ CART-A CLEANING AGENT

# DIADEZ CART -A

## cartridges

With DIADEZ CART-A cartridges, regular dialysis machine maintenance and disinfection assure proper dialysis machine function and requires less intensive machine operation from clinical staff.

- Less space for storage and less heavy lifting for staff compared to liquid disinfectant
- Minimizes risk of personal injury and material damage when handling chemical agents
- Easy to streamline disinfection combining heat disinfection with decalcification and cleaning
- DIADEZ CART-A cartridge removes organic deposits, fats and proteins in combination with the heat disinfection program of the machine. It is recommended at least after 22 treatments or at least once every seventh treatment day
- Recyclable polypropylene container for less environmental impact compared to the liquid disinfectant

Dialysis machine maintenance and disinfection becomes a convenient and simple operation.





## **DIADEZ CART-A CARTRIDGE SPECIFICATIONS**

The **DIADEZ CART-A** cartridge contains 13 grams of anhydrous sodium carbonate powder. The **DIADEZ CART-A** cartridge is intended for preparation of a sodium carbonate solution used for removing organic deposits, fats and proteins from the dialysis machine's fluid circuit. The **DIADEZ CART-A** cartridge must be used in combination with the heat disinfection program of the corresponding dialysis machine.

TECHNICAL DATA	DIADEZ CART-A
Action in combination with the heat disinfection program	Cleaning of the fluid path from organic deposits, fats, and proteins
Cartridge content	13 grams
Total cartridge weight	Approx. 32 grams
Cartridge and label material	Recyclable polypropylene (PP)
Recommended usage	Refer to operator's manual for corresponding dialysis machine
pH (prepared solution)	Approx. 11
Residual test needed	No
Biodegradable content	Yes
Storage conditions	Store below +30°C
Shelf life	24 months from date of manufacture
Packaging	Box of 80 cartridges

PRODUCT	PACKAGING	UNIT PER PACKAGE	PACKAGE SIZE (cm)	PACKAGE WEIGHT (kg)	PALLET SIZE (cm)	PALLET WEIGHT (kg)	UNIT PER PALLET	BOX PER PALLET
Diadez Cart A	Box	80	25*40*18	3	80*120*156	245	5760	72



MEDICA PONTEM



MEDICA PONTEM



## DIADÉZ CART-32 CARTRIDGE SPECIFICATIONS



MEDICA PONTEM

## DIADEZ CART-32 CARTRIDGE SPECIFICATIONS

### DIADEZ CART 32 g

#### Disinfectant for GAMBRO Hemodialysis Machines

The DIADEZ Cart cartridge contains 32 grams of citric acid powder. The DIADEZ CART cartridge is intended for preparation of a citric acid solution used for removing calcium and magnesium precipitation from the dialysis machine's fluid circuit. The DIADEZ CART cartridge must be used in combination with the heat disinfection program of the corresponding dialysis machine.

TECHNICAL DATA	DIADEZ CART
Action in combination with the heat disinfection program	Decalcification of the fluid path
Cartridge content	32 grams
Total cartridge weight	Approx. 50 grams
Cartridge and label material	Recyclable polypropylene (PP)
Recommended usage	Refer to operator's manual for corresponding dialysis machine
pH (prepared solution)	Approx. 2
Residual test needed	No
Biodegradable content	Yes
Storage conditions	Store below +30°C
Shelf life	36 months from date of manufacturing date
Packaging	Box of 33 cartridges
Box dimensions	279 x 247 x 81 mm
Total box weight	1800 grams

PRODUCT	PACKAGING	UNIT PER PACKAGE	PACKAGE SIZE (cm)	PACKAGE WEIGHT (kg)	PALLET SIZE (cm)	PALLET WEIGHT (kg)	UNIT PER PALLET	BOX PER PALLET
Diadez Cart 32	Box	80	25*40*18	4,5	80*120*156	350	5760	72



MEDICA PONTEM